

OF-1 Technical Data & Information Sheet

1- Origin

OF-1 is extracted from a seaweed called Funori.

It is 100% made in Japan.

In particular, we use the ones from the Kyusyu region in Japan, which are said to have strong glue.

2- Composition

It is a polysaccharide extracted from seaweed Funori.

3- How to Use

(1) Weight only the amount of OF-1 you want to use.

(2) Dissolve in water.

The standard formula is to dissolve 1g of OF-1 in 100g of water.

After that, adjust the concentration to your liking.

(3) Raise the temperature while stirring.

Please do not bring it to a boil.

Keep the temperature up to 80°C.

(4) After dissolve, cool it naturally before use.

***Attention** Freeze-dried products may clump if they are melted suddenly from hot water or if the stirring time is short. In that case, stir well and check that there are no lump before using.

4- Package

OF-1 is packed in an aluminum zipper bag.

It is 1g per bag.

※The package is currently only 1g.

5- Treatment & Storage

Avoid direct sunlight and store in a cool and dark place.

There is no problem at room temperature.

And, If unopened, there is no problem with quality for several years(Result of our Test)

Owaki Manzo LTD.,CO

info@funori.com